

LIST OF TABLES

TABLE NO.	TITLE	PAGE NO.
1.1	Dielectric properties of the solvents used commonly for Microwave extraction	23
2.1	Optimisation of extractant volume with respect to weight of solid sample to achieve the best recovery. Recoveries of PCP for soil sample spiked with standards, at 0.5 $\mu\text{g/g}$ and leather sample (originally having PCP at 1.3 $\mu\text{g/g}$). Detection by GC-ECD. MAE conditions: 20% MWP for 15min duration)	69
2.2	LOD of chlorophenols found by spiking the standards, MAE with 0.01% carbonate solution and ENVI-18 cleanup (by HPLC-UV); and derivatised before ENVI-18 cleanup for GC-ECD detection. Recovery (%) of chlorophenols found by spiking standards at 1 $\mu\text{g/g}$ (each analyte) to solid samples by aqueous MAE and ENVI-18 follow-up (% RSD 5 replicates)	73
2.3	Recovery values of chlorophenols of different spike levels, from solid samples extracted with 0.01% aqueous carbonate by MW, acetyl derivatised followed by ENVI-18 cleanup. (Detection by GC-ECD % RSD from 3 replicates)	74

TABLE NO.	TITLE	PAGE NO.
2.4	Detection (by GC-ECD) of PCP in real samples by different extraction techniques	75
3.1	LOD and % Recoveries (ng/g) for spiked analytes of soil samples	110
3.2	LOD and % Recoveries (ng/ml) for spiked analytes of water samples	110
3.3	Comparison of the recovery values obtained for the spiked samples analyzed by various methods. %RSD (n=3) is given in parenthesis	114
4.1	LOD and % Recoveries (ng/g) for spiked analytes from solid samples	128
4.2	LOD and % Recoveries (ng/ml) for spiked analytes of water samples	132
4.3	Detection (by GC-ECD) of PCP in real samples by different extraction techniques	142
5.1	Recoveries for chlorophenols by MASD with different on-line enrichments followed	166
5.2	Recoveries for chloro organic compounds by MASD with different on-line enrichments followed	166
5.3	Comparison of recoveries by different SPE sorbents used for on-line enrichments	171
5.4	Results obtained for real samples and spiked samples by MASD with different on-line enrichments	173
5.5	LOD (ng/g) for spiked analytes in solid (soil) samples (calculated as 5 times the base line noise signal)	174