

Appendix V

Composition of AB medium

20X AB buffer	K ₂ HPO ₄ NaH ₂ PO ₄ Water Sterilization	60g 20g 1000ml Autoclaving
20X AB salt	NH ₄ Cl MgSO ₄ . 7H ₂ O KCl CaCl ₂ FeSO ₄ . 7H ₂ O Water Sterilization	20g 6g 3g 0.2g 0.04g 1000ml Autoclaving
	Glucose Water	5g 900ml
AB media	AB salt AB buffer Glucose	50ml 50ml 900ml
All should be added after autoclaving under		